

YMC Chiral Applications

S.No	Name of the Chiral Molecule	Name of the Column & Dimension
1	Abacavir USP	YMC Chiral art Amylose C Neo (4.6x250mm, 5µm)
2	Acalabrutinib	Sumichiral OA-5000L (4.6x150mm, 5µm)
3	Adrenaline (Epinephrine)	SUMI Chiral OA-6100(150X4.6mm) 5µm
4	Alanamide	YMC Chiral art Cellulose SZ (4.6x250mm, 5µm)
5	Alanine	YMC Chiral art Cellulose SB (4.6x250mm, 5µm)
6	Albaconazol	YMC Chiral art Cellulose SC (4.6x250mm, 5µm)
7	all-trans-Astaxanthin	YMC Chiral art Amylose SA (4.6x250mm, 5µm)
8	Amlodipine	YMC Chiral art Cellulose SC (4.6x250mm, 5µm)
9	Apixaban	YMC Chiral art Cellulose SZ (4.6x150mm, 3µm)
10	Aprepitant	YMC Chiral art Cellulose SZ (4.6x250mm, 5µm)
11	Aprepitant	YMC Chiral art Amylose C Neo (4.6x150mm, 5µm)
12	Argatroban	YMC Chiral art Amylose C (4.6x250mm, 5µm)
13	Atorvastatin EP	YMC Chiral art Cellulose C (4.6x250mm, 5µm)
14	Atorvastatin USP	YMC Chiral art Cellulose C (4.6x250mm, 5µm)
15	Atropine	YMC Chiral art Amylose C (4.6x250mm, 5µm)
16	Atropine	YMC Chiral art Cellulose C (4.6x250mm, 5µm)
17	Atropine	YMC Chiral art Amylose SA (4.6x250mm, 5µm)
18	Avacopan	YMC Chiral art Amylose C Neo (4.6x250mm, 5µm)
19	Bedaquiline Fumarate	YMC Chiral art Cellulose SC (3x250mm, 3µm)
20	Benidipine	YMC Chiral art Cellulose SZ (4.6x250mm, 5µm)
21	Benzoin	YMC Chiral art Cellulose C (4.6x150mm, 3µm)
22	Benzoin	YMC Chiral art Cellulose SJ (4.6x250mm, 5µm)
23	Bepotastine	YMC Chiral art Amylose SA (4.6x250mm, 5µm)
24	Brinzolamide	YMC Chiral art Cellulose SB (4.6x250mm, 5µm)
25	Brivaracetam	YMC Chiral art Cellulose C (4.6x250mm, 5µm)
26	Brivaracetam	YMC Chiral art Amylose SA (3x250mm, 3µm)
27	Bromuconazole	YMC Chiral art Cellulose SZ (4.6x150mm, 3µm)
28	Cabotegravir sodium	YMC Chiral art Cellulose C (3x150mm, 3µm)
29	Canagliflozin	YMC Chiral art Cellulose SZ (4.6x250mm, 5µm)
30	Carprofen	YMC Chiral art Cellulose SJ (4.6x250mm, 5µm)
31	Cenobamate	YMC Chiral art Cellulose SB (4.6x150mm, 3µm)
32	Cetirizine (Levocetirizine)	YMC Chiral art Amylose SA (3x250mm, 3µm)
33	Cetirizine (Levocetirizine)	YMC Chiral art Cellulose SC (4.6x250mm, 5µm)
34	Chlormezanone	YMC Chiral art Amylose SA (4.6x250mm, 5µm)
35	Chloroquine	YMC Chiral art Cellulose SC (4.6x250mm, 3µm)
36	Chrysanthemic acid	YMC Chiral art Amylose C Neo (4.6x150mm, 3µm)
37	Citalopram	YMC Chiral art Amylose C Neo (4.6x250mm, 5µm)
38	Citalopram hydrobromide	YMC Chiral art Cellulose SB (4.6x250mm, 5µm)
39	Clopidogral	YMC Chiral art Amylose SA (3x250mm, 3µm)
40	Clopidogrel USP	YMC Chiral art Cellulose SJ (4.6x250mm, 5µm)
41	Colchicine	YMC Chiral art Amylose C (4.6x250mm, 5µm)
42	Coriandrol Licareol	YMC Chiral art Cellulose C (4.6x250mm, 5µm)
43	Cyproconazole	YMC Chiral art Cellulose SB (4.6x250mm, 5µm)
44	Darolutamide	YMC Chiral art Cellulose SC (2x150mm, 3µm)
45	Darunavir	YMC Chiral art Cellulose C (4.6x250mm, 5µm)
46	Desglymidodrine	YMC Chiral art Cellulose SB (4.6x250mm, 5µm)

47	Dexketoprofen	YMC Chiral art Amylose C (4.6x250mm, 5um)
48	Dexmedetomidine	YMC Chiral art Cellulose C (4.6x250mm, 5um)
49	Diclobutazole	YMC Chiral art Cellulose SZ (4.6x150mm, 5um)
50	Diniconazole	YMC Chiral art Cellulose SB (4.6x250mm, 5um)
51	DL-Aminoglutethimide	YMC Chiral art Cellulose SJ (4.6x250mm, 5um)
52	DL-Mandelic acid	YMC Chiral art Cellulose SZ (4.6x250mm, 5um)
53	DL-Phenylethyl amine	YMC Chiral art Cellulose SB (3x150mm, 3um)
54	Dolutegravir	YMC Chiral art Cellulose SB (4.6x250mm, 5um)
55	Donepezil	YMC Chiral art Cellulose SZ (4.6x250mm, 5um)
56	Duloxetine Hydrochloride	YMC Chiral art Cellulose SJ (4.6x250mm, 5um)
57	Duloxetine Hydrochloride USP	YMC Chiral art Amylose SA (4.6x250mm, 5um)
58	Econazole	YMC Chiral art Cellulose SA (4.6x250mm, 5um)
59	Edoxaban	YMC Chiral art Amylose C (4.6x250mm, 5um)
60	Elacestrant	YMC Chiral art Cellulose SZ (4.6x250mm, 5um)
61	Eletriptan hydrobromide	Sumichiral OA-6100 (4.6x150mm, 5um)
62	Empagliflozin	YMC Chiral art Cellulose SZ (4.6x150mm, 5um)
63	Empagliflozin	YMC Chiral art Cellulose SC (3x250mm, 3um)
64	Enilconazole	YMC Chiral art Cellulose C (4.6x250mm, 5um)
65	Epoconazole	YMC Chiral art Amylose SA (4.6x250mm, 5um)
66	Equol	YMC Chiral art Amylose SA (4.6x250mm, 3um)
67	Escitalopram	YMC Chiral art Cellulose SJ (4.6x250mm, 5um)
68	Esomeprazole	YMC Chiral art Amylose C Neo (4.6x250mm, 5um)
69	Famprofazone	YMC Chiral art Cellulose SB (4.6x250mm, 5um)
70	Fenoprofen	YMC Chiral art Cellulose C (4.6x250mm, 5um)
71	Fexofenadine	YMC Chiral art Amylose C Neo (4.6x250mm, 5um)
72	Florfenicol	YMC Chiral art Amylose SJ (4.6x250mm, 5um)
73	Flurbiprofen	YMC Chiral art Cellulose SJ (4.6x250mm, 5um)
74	Formoterol	YMC Chiral art Cellulose SC (4.6x250mm, 5um)
75	Fulvestrant	YMC Chiral art Cellulose SJ (4.6x250mm, 5um)
76	Galantamine	YMC Chiral art Cellulose SJ (4.6x250mm, 5um)
77	Glutamine	YMC Chiral art Amylose C Neo (4.6x250mm, 5um)
78	Hexaconazole	YMC Chiral art Amylose SA (4.6x250mm, 5um)
79	Hexaconazole	YMC Chiral art Cellulose SC (3x250mm, 3um)
80	Hydroxychloroquine	YMC Chiral art Cellulose SB (4.6x150mm, 3um)
81	Ibrutinib	YMC Chiral art Cellulose SA (4.6x250mm, 5um)
82	Ibuprofen	YMC Chiral art Amylose C Neo (4.6x250mm, 5um)
83	Imiprothrin	YMC Chiral art Amylose C Neo (4.6x250mm, 5um)
84	Indacaterol	YMC Chiral art Cellulose C (4.6x250mm, 10um)
85	Indapamide	YMC Chiral art Cellulose C (4.6x250mm, 5um)
86	Isavuconazole	YMC Chiral art Cellulose SB (4.6x150mm, 3um)
87	Isavuconazole	YMC Chiral art Amylose C Neo (4.6x250mm, 5um)
88	Isavuconazonium sulfate	YMC Chiral art Amylose C Neo (4.6x250mm, 5um)
89	Isoxsuprine	YMC Chiral art Amylose SA (4.6x250mm, 5um)
90	Ketamine	YMC Chiral art Amylose SA (4.6x250mm, 5um)
91	Ketoprofen	YMC Chiral art Cellulose C (4.6x150mm, 3um)
92	Ketorolac	YMC Chiral art Cellulose C (4.6x250mm, 5um)
93	Ketorolac tromethamine	YMC Chiral art Cellulose SB (4.6x250mm, 5um)
94	Lacosamide USP	YMC Chiral art Cellulose SB (4.6x250mm, 5um)
95	Lansoprazole	YMC Chiral art Amylose C Neo (4.6x150mm, 3um)
96	Lanzoprazole	YMC Chiral art Amylose C (4.6x250mm, 5um)
97	Latnoprost	YMC Chiral art Cellulose C (4.6x250mm, 5um)
98	Levetiracetam USP method	YMC Chiral art Cellulose SB (4.6x250mm, 5um)
99	Levetiracetam-EP Method	YMC Chiral art Cellulose SB (4.6x250mm, 5um)
100	Levocetirizine	YMC Chiral art Cellulose SB (4.6x250mm, 5um)

101	Levocetirizine	YMC Chiral art Cellulose C (4.6x250mm, 5um)
102	Levonorgestrel	YMC Chiral art Cellulose SC (4.6x250mm, 5um)
103	Levosulpiride	YMC Chiral art Cellulose SC (4.6x250mm, 5um)
104	Levosulpiride	YMC Chiral art Cellulose SC (4.6x250mm, 5um)
105	Linagliptin	YMC Chiral art Amylose C (4.6x250mm, 5um)
106	Linezolid	YMC Chiral art Cellulose SC (4.6x250mm, 3um)
107	Luliconazole	YMC Chiral art Cellulose SC (4.6x250mm, 5um)
108	Luliconazole Cream Assay	YMC Chiral art Cellulose SC (4.6x250mm, 5um)
109	Luliconazole Drug Substance Assay	YMC Chiral art Amylose SA (3x250mm, 3um)
110	Lurasidone	YMC Chiral art Amylose SA (4.6x250mm, 5um)
111	Metalaxyl	YMC Chiral art Cellulose SZ (4.6x150mm, 5um)
112	Metconazole	YMC Chiral art Cellulose SB (3x150mm, 3um)
113	Metoprolol	YMC Chiral art Cellulose SJ (4.6x250mm, 5um)
114	Metoprolol	YMC Chiral art Cellulose SJ (4.6x250mm, 5um)
115	Metoprolol	YMC Chiral art Cellulose SJ (4.6x250mm, 5um)
116	Metoprolol tartrate	YMC Chiral art Cellulose SC (4.6x250mm, 5um)
117	Mirogabalin	YMC Chiral art Cellulose SC (4.6x250mm, 5um)
118	Naproxen	YMC Chiral art Amylose SA (4.6x150mm, 3um)
119	Naringenin	YMC Chiral art Amylose SA (4.6x250mm, 3um)
120	Nelarabine	YMC Chiral art Amylose SA (4.6x250mm, 3um)
121	Nelarabine	YMC Chiral art Cellulose-SC (4.6x250mm, 3um)
122	Netarsudil	Sumichiral OA-5000L (4.6x150mm, 5um)
123	Netarsudil	YMC Chiral art Cellulose C (3x150mm, 3um)
124	Nimodipine	YMC Chiral art Cellulose C (4.6x250mm, 5um)
125	Noradrenaline (Norepinephrine)	YMC Chiral art Cellulose C (4.6x250mm, 5um)
126	Ofloxacin (Levofloxacin)	YMC Chiral art Cellulose SZ (4.6x250mm, 5um)
127	Ofloxacin (Levofloxacin)	YMC Chiral art Cellulose SC (4.6x250mm, 5um)
128	Omeprazole	YMC Chiral art Cellulose SC (4.6x250mm, 5um)
129	Oxybutynin	YMC Chiral art Cellulose SC (4.6x250mm, 5um)
130	Paclobutrazol	YMC Chiral art Cellulose SB (4.6x250mm, 5um)
131	Paclobutrazol	YMC Chiral art Amylose SA (4.6x250mm, 5um)
132	Phenoxybenzamine	YMC Chiral art Amylose SA (4.6x250mm, 3um)
133	Phenoxybenzamine hydrochlorid	YMC Chiral art Cellulose SB (4.6x250mm, 5um)
134	Phenoxybenzamine	YMC Chiral art Cellulose-SC (4.6x250mm, 5um)
135	Phenylephrine	YMC Chiral art Amylose C (4.6x250mm, 5um)
136	Phenylethylamine	YMC Chiral art Amylose C (4.6x250mm, 5um)
137	Phytanodione	YMC Chiral art Amylose C Neo (4.6x250mm, 5um)
138	Pidotimod	YMC Chiral art Amylose SA (4.6x250mm, 5um)
139	Pindolol	YMC Chiral art Amylose C Neo (4.6x250mm, 5um)
140	Pindolol	YMC Chiral art Amylose SA (4.6x250mm, 3um)
141	Pindolol	YMC Chiral art Amylose SA (3x250mm, 3um)
142	Pioglitazone	YMC Chiral art Amylose C (4.6x250mm, 5um)
143	Ponesimod	YMC Chiral art Amylose C Neo (4.6x250mm, 5um)
144	Ponesimod	YMC Chiral art Cellulose SC (3x250mm, 3um)
145	Proline	YMC Chiral art Amylose C (3x150mm, 3um)
146	Propiconazole	YMC Chiral art Cellulose SJ (4.6x250mm, 5um)
147	Propranolol	YMC Chiral art Cellulose SC (4.6x250mm, 3um)
148	Propranolol	YMC Chiral art Cellulose SB (4.6x250mm, 5um)
149	Pyroglutamic Acid	YMC Chiral art Cellulose SZ (4.6x250mm, 5um)
150	Rabeprazole	YMC Chiral art Cellulose C (3x150mm, 3um)
151	Resveratrol	YMC Chiral art Amylose C (3x150mm, 3um)
152	Rimegepant	YMC Chiral art Cellulose SJ (4.6x250mm, 5um)
153	Rivaroxaban	YMC Chiral art Cellulose SJ (4.6x250mm, 5um)
154	Rivaroxaban	YMC Chiral art Cellulose SB (3x150mm, 3um)

155	Ropivacaine	YMC Chiral art Amylose SA (4.6x250mm, 5um)
156	Rosiglitazone	YMC Chiral art Cellulose SC (4.6x250mm, 5um)
157	Rosuvastatin	YMC Chiral art Amylose SA (4.6x250mm, 5um)
158	Ruxolitinib	YMC Chiral art Cellulose C (4.6x250mm, 5um)
159	Sertraline hydrochloride EP	YMC Chiral art Cellulose C (4.6x250mm, 5um)
160	Sertraline hydrochloride SH(enantiomer)	YMC Chiral art Amylose C (4.6x250mm, 5um)
161	Sitagliptin	YMC Chiral art Cellulose SB (4.6x250mm, 5um)
162	Sitagliptin	YMC Chiral art Cellulose SJ (3x50mm, 3um)
163	Sitagliptin USP	YMC Chiral art Cellulose SC (4.6x250mm, 5um)
164	Sofpironium Bromide	YMC Chiral art Cellulose SZ (4.6x250mm, 5um)
165	Sotorasib	YMC Chiral art Cellulose SB (4.6x250mm, 5um)
166	Spiroxamine	YMC Chiral art Cellulose SC (2x150mm, 3um)
167	Suvorexant	YMC Chiral art Cellulose SC (4.6x250mm, 5um)
168	Tadalafil	YMC Chiral art Amylose C Neo (4.6x250mm, 5um) Dought
169	Tapentadol	YMC Chiral art Cellulose SC (4.6x250mm, 5um)
170	Tebuconazol	YMC Chiral art Cellulose SC (4.6x250mm, 5um)
171	Tetrahydrocannabinol (THC)	YMC Chiral art Cellulose SZ (4.6x250mm, 5um)
172	Thalidomide	YMC Chiral art Amylose C Neo (4.6x250mm, 5um)
173	Tiagabine	YMC Chiral art Cellulose SC (4.6x250mm, 5um)
174	Tofacifinib	YMC Chiral art Cellulose SZ (4.6x250mm, 5um)
175	Tofacitinib	YMC Chiral art Cellulose SZ (4.6x250mm, 5um)
176	Trelagliptin	YMC Chiral art Cellulose SC (4.6x250mm, 5um)
177	Triadimefon	YMC Chiral art Cellulose SZ (4.6x250mm, 5um)
178	Triadimenol	YMC Chiral art Amylose SA (4.6x250mm, 5um)
179	Trimebutine	YMC Chiral art Cellulose SC (4.6x250mm, 5um)
180	Trimebutine maleate	YMC Chiral art Amylose SA (4.6x250mm, 5um)
181	Triticonazole	YMC Chiral art Cellulose SC (4.6x250mm, 5um)
182	Troger's base	YMC Chiral art Cellulose SC (4.6x250mm, 5um)
183	Ubrogepant	YMC Chiral art Cellulose SZ (4.6x250mm, 5um)
184	Udenafil	YMC Chiral art Cellulose SC (4.6x250mm, 5um)
185	Valsartan (USP)	YMC Chiral art Amylose C Neo (4.6x250mm, 5um)
186	Valsartan Diclofenac sodium (I.S.)	YMC Chiral art Amylose SA (4.6x250mm, 5um)
187	Verapamil	YMC Chiral art Cellulose SJ (4.6x250mm, 5um)
188	Vonoprazan	YMC Chiral art Amylose SA (4.6x250mm, 5um)
189	Warfarin	YMC Chiral art Cellulose SB (4.6x250mm, 5um)
190	Warfarin	YMC Chiral art Cellulose SC (4.6x250mm, 5um)
191	Zanubrutinib	YMC Chiral art Cellulose SC (4.6x250mm, 5um)
192	Zanubrutinib	YMC Chiral art Cellulose SB (4.6x250mm, 5um)
193	Zanubrutinib	YMC Chiral art Amylose C Neo (4.6x250mm, 5um)
194	Zopiclone	YMC Chiral art Cellulose SC (4.6x250mm, 5um)



YMC Application LAB

Second Floor, Plot No- 78/B3,
Phase VI, IDA Jeedimetla,
Quthbullapur Mandal,
Medchal-Malkajgiri District,
Telangana-500055



HEAD OFFICE

Ground Floor-01, C-125 A,
Block C, Sector-2,Phase-1, Noida,
Gautambuddha Nagar,(U.P.)-201301.
Tel: 0120-4369405
E-mail : sales@ymcindia.com
Website : www.ymcindia.com



+91 8826239666



+91 8826239666



04029303386 / 236



applicationlabsport@ymcindia.com